



	Height 1	ID Barrel 2	Fin 3	OD Barrel 4	OD Septum 5	Top 6	Weight 7
20mm	16.50 +/- .15	9.50 +/- .15	0.60 +/- .1	11.35 +/- .15	8.10 +/- .15	17.10 +/- .15	0.80 +/- .1
24mm	17.30 +/- .15	14.1 +/- .15	0.70 +/- .1	16.50 +/- .15	8.10 +/- .15	21.0 +/- .15	1.40 +/- .1
28mm	17.70 +/- .15	15.65 +/- .15	0.70 +/- .1	18.70 +/- .15	8.40 +/- .15	24.6 +/- .15	1.85 +/- .1
	Min	Mid	Max				
20mm	12.6	13.9	15.3				
24mm	17.9	18.7	19.5				
28mm	20.2	20.8	21.4				

These dimensions are guidelines only.

SealSafe® adapters should be tested for use with specific containers, liquids and oral syringe tips to confirm fit and compatibility and to ensure optimal performance.

The Andwin Corp. shall not be liable for any consequential damages, loss of liquids, loss of or damage to property or claims of third parties arising out of or in connection with use of SealSafe® adapters.